


2-CONnect®

Angle compensation for divergent implants available in different gingiva heights and angled version



BENEFITS:

- self - centering conical connection
- maximum flexibility of divergent implant positions
- stressless fitting by optional adhesive connection (female part)
- force application caused by tube shape and integrated screw (Straight Abutment) or mounted screw (Angled Abutment)
- female part compatible between the different implant systems and diameters



APPLICATION:

- bars
- bridges
- non-engaging structures



TECHNICAL DETAILS:

- Titanium Grade 5 ELI

COPING (FEMALE PART) ON 2-CONNECT®

- *Female = glued



ABUTMENT (MALE PART) ON 2-CONNECT®

- *Male = screwed

