


COL-SERIES

compatible to:

Altatec[®]

Conelog[®]



A: over surface

B: over edge

PLATFORMS:



3,3 mm



3,8 / 4,3 mm



5,0 mm



TITANIUM BASE

incl. screws

GH 0,8 mm

COL 80330008

GH 2,0 mm

COL 80330020



SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

COL 9.S3D3.300



NT-MILL LAB ANALOG

COL 5.La3.300



IMPLANT PICK-UP

for open-tray with integrated screw

COL TR-00013.3



COLOURED LAB SCREW

COL 623343



ABUTMENT SCREW



Overview of all screwdrivers on pages 344 and 345

COL 603343



DIM ANALOG®

COL 5.DIM.330



SLEEVE (belongs to DIM ANALOG)

plaster model

DIM 5.XSL.A10





<p>NT-PREFORM® incl. coloured milling screw</p>	<table border="1"> <tr> <td data-bbox="798 403 933 448">GH 0,5 mm</td> <td data-bbox="933 403 1220 448">COL 10.Mi3.300.05</td> </tr> <tr> <td data-bbox="798 448 933 504">GH 1,3 mm</td> <td data-bbox="933 448 1220 504">COL 10.Mi3.300</td> </tr> </table>	GH 0,5 mm	COL 10.Mi3.300.05	GH 1,3 mm	COL 10.Mi3.300	
GH 0,5 mm	COL 10.Mi3.300.05					
GH 1,3 mm	COL 10.Mi3.300					
<p>SCAN BODY 3D-GUIDE® for intra + extra oral scanning</p>	<p>COL 9.S3D3.300</p>					
<p>NT-MILL LAB ANALOG for NT-Preform®</p>	<p>COL 5.La3.300</p>	<p>groove</p>				
<p>COLOURED MILLING SCREW for NT-Preform®</p>	<p>COL 6.Sc3.343</p>					
<p>ABUTMENT SCREW</p>	<p>COL 603343</p>					
<p>Overview of all screwdrivers on pages 344 and 345</p>						
<p>DIM ANALOG®</p>	<p>COL 5.DIM.330</p>					
<p>SLEEVE (belongs to DIM ANALOG) plaster model</p>	<p>DIM 5.XSL.A10</p>					
<p>ABUTMENT GRIP for NT-Holder fixture</p>	<p>W 12.100.000</p>					
<p>CUSTOMIZED ABUTMENT milled by NT-iDent</p>	<p>Not available in the USA TZ15.TiO.Si1</p>					

COL
SERIES



2-CONNECT®-ABUTMENT SET (STRAIGHT)

incl. 2-CONnect®-Coping and Coping-screw

GH 0,5 mm COL-2CON 8033001S

GH 2,5 mm COL-2CON 8033002S



2-CONNECT®-ABUTMENT

Male Part

GH 0,5 mm COL-2CON 8033001M

GH 2,5 mm COL-2CON 8033002M



2-CONNECT®-SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

KS 61 = 1,8 mm

2CON 9.S3D1.800



2-CONNECT®-COPING

Female Part

2CON F-U



COPING-SCREW

Recommended tightening torques: 15 Ncm

KS 61



2-CONNECT®-DIM ANALOG®

KS 61 = 1,8 mm

2CON 5.DIM.180



SLEEVE (belongs to DIM ANALOG)

plaster model

DIM 5.XSL.A30



2-CONNECT®-IMPLANT PICK-UP

integrated screw

KS 61 = 1,8 mm

2CON TR-00021.8



2-CONNECT®-LAB ANALOG

KS 61 = 1,8 mm

2CON 5.La1.800





<p>2-CONNECT®-ABUTMENT SET (ANGLED) incl. 2-CONnect®-Coping, mounting aid Coping-screw and Male Part Screw, angled 20°, Gingiva Height 2,5 mm</p>	<p>A: over surface B: over edge</p>	<p>COL-2CON 80330C3S COL-2CON 80330A3S</p>	
<p>2-CONNECT®-ABUTMENT angled 20°, Gingiva Height 2,5 mm Male Part incl. screw, mounting aid</p>	<p>A: over surface B: over edge</p>	<p>COL-2CON 80330C3M COL-2CON 80330A3M</p>	
<p>2-CONNECT®-SCAN BODY 3D-GUIDE® for intra + extra oral scanning KS 62 = 1,6 mm</p>		<p>2CON 9.S3D1.600</p>	
<p>2-CONNECT®-COPING Female Part</p>		<p>2CON F-U</p>	
<p>COPING-SCREW Recommended tightening torques: 15 Ncm</p>		<p>KS 62</p>	
<p>MALE PART SCREW for angled 20° 2-CONnect® Abutment Recommended tightening torques: 20 Ncm</p>		<p>COL 6.Ms3.343</p>	
<p>2-CONNECT®-DIM ANALOG® KS 62 = 1,6 mm</p>		<p>2CON 5.DIM.160</p>	
<p>SLEEVE (belongs to DIM ANALOG) plaster model</p>		<p>DIM 5.XSL.A30</p>	
<p>2-CONNECT®-IMPLANT PICK-UP integrated screw KS 62 = 1,6 mm</p>		<p>2CON TR-00031.6</p>	
<p>2-CONNECT®-LAB ANALOG KS 62 = 1,6 mm</p>		<p>2CON 5.La1.600</p>	



TITANIUM BASE

incl. screws

GH 0,8 mm	COL 80384308
GH 2,0 mm	COL 80384320



SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

COL 9.S3D3.843



NT-MILL LAB ANALOG

COL 5.La3.843



IMPLANT PICK-UP

for open-tray with integrated screw

COL TR-00023.8



COLOURED LAB SCREW

COL 623343



ABUTMENT SCREW



Overview of all screwdrivers on pages 344 and 345

COL 603343



DIM ANALOG®

COL 5.DIM.384



SLEEVE (belongs to DIM ANALOG)
plaster model

DIM 5.XSL.A20





<p>NT-PREFORM® incl. coloured milling screw</p>	<p>GH 0,5 mm COL 10.Mi3.843.05 GH 1,3 mm COL 10.Mi3.843</p>	
<p>SCAN BODY 3D-GUIDE® for intra + extra oral scanning</p>	<p>COL 9.S3D3.843</p>	
<p>NT-MILL LAB ANALOG for NT-Preform®</p>	<p>COL 5.La3.843</p>	
<p>COLOURED MILLING SCREW for NT-Preform®</p>	<p>COL 6.Sc3.343</p>	
<p>ABUTMENT SCREW</p>	<p>COL 603343</p>	
<p>Overview of all screwdrivers on pages 344 and 345</p>		
<p>DIM ANALOG®</p>	<p>COL 5.DIM.384</p>	
<p>SLEEVE (belongs to DIM ANALOG) plaster model</p>	<p>DIM 5.XSL.A20</p>	
<p>ABUTMENT GRIP for NT-Holder fixture</p>	<p>W 12.200.000</p>	
<p>CUSTOMIZED ABUTMENT milled by NT-iDent</p>	<p>Not available in the USA TZ15.TiO.Si1</p>	

COL
SERIES



2-CONNECT®-ABUTMENT SET (STRAIGHT)

incl. 2-CONnect®-Coping and Coping-screw

GH 0,5 mm COL-2CON 8038431S

GH 2,5 mm COL-2CON 8038432S



2-CONNECT®-ABUTMENT

Male Part

GH 0,5 mm COL-2CON 8038431M

GH 2,5 mm COL-2CON 8038432M



2-CONNECT®-SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

KS 61 = 1,8 mm

2CON 9.S3D1.800



2-CONNECT®-COPING

Female Part

2CON F-U



COPING-SCREW

Recommended tightening torques: 25 Ncm

KS 61



2-CONNECT®-DIM ANALOG®

KS 61 = 1,8 mm

2CON 5.DIM.180



SLEEVE (belongs to DIM ANALOG)

plaster model

DIM 5.XSL.A30



2-CONNECT®-IMPLANT PICK-UP

integrated screw

KS 61 = 1,8 mm

2CON TR-00021.8



2-CONNECT®-LAB ANALOG

KS 61 = 1,8 mm

2CON 5.La1.800






2-CONNECT®-ABUTMENT SET (ANGLED)

incl. 2-CONnect®-Coping, mounting aid
Coping-screw and Male Part Screw, angled 20°,
Gingiva Height 2,5 mm

A: over surface	COL-2CON 80384C3S
B: over edge	COL-2CON 80384A3S



2-CONNECT®-ABUTMENT

angled 20°, Gingiva Height 2,5 mm
Male Part incl. screw,
mounting aid

A: over surface	COL-2CON 80384C3M
B: over edge	COL-2CON 80384A3M



2-CONNECT®-SCAN BODY 3D-GUIDE®

for intra + extra oral scanning
KS 62 = 1,6 mm

2CON 9.S3D1.600



2-CONNECT®-COPING

Female Part

2CON F-U



COPING-SCREW

Recommended tightening torques: 15 Ncm

KS 62



MALE PART SCREW

for angled 20° 2-CONnect® Abutment
Recommended tightening torques: 20 Ncm

COL 6.Ms3.343



2-CONNECT®-DIM ANALOG®

KS 62 = 1,6 mm

2CON 5.DIM.160



SLEEVE (belongs to DIM ANALOG)

plaster model

DIM 5.XSL.A30



2-CONNECT®-IMPLANT PICK-UP

integrated screw
KS 62 = 1,6 mm

2CON TR-00031.6



2-CONNECT®-LAB ANALOG

KS 62 = 1,6 mm

2CON 5.La1.600





TITANIUM BASE

incl. screws

GH 0,8 mm

COL 80500008

GH 2,0 mm

COL 80500020



SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

COL 9.S3D5.000



NT-MILL LAB ANALOG

COL 5.La5.000



COLOURED LAB SCREW

COL 625000



ABUTMENT SCREW



Overview of all screwdrivers on pages 344 and 345

COL 605000



DIM ANALOG®

COL 5.DIM.500



SLEEVE (belongs to DIM ANALOG)

plaster model

DIM 5.XSL.A30



<p>NT-PREFORM® incl. coloured milling screw</p>	<p>GH 0,5 mm COL 10.Mi5.000.05 GH 1,3 mm COL 10.Mi5.000</p>	
<p>SCAN BODY 3D-GUIDE® for intra + extra oral scanning</p>	<p>COL 9.S3D5.000</p>	
<p>NT-MILL LAB ANALOG for NT-Preform®</p>	<p>COL 5.La5.000</p>	 groove
<p>COLOURED MILLING SCREW for NT-Preform®</p>	<p>COL 6.Sc5.000</p>	
<p>ABUTMENT SCREW</p>  <p>Overview of all screwdrivers on pages 344 and 345</p>	<p>COL 605000</p>	
<p>DIM ANALOG®</p>	<p>COL 5.DIM.500</p>	
<p>SLEEVE (belongs to DIM ANALOG) plaster model</p>	<p>DIM 5.XSL.A30</p>	
<p>ABUTMENT GRIP for NT-Holder fixture</p>	<p>W 12.300.000</p>	
<p>CUSTOMIZED ABUTMENT milled by NT-iDent</p> <p>Not available in the USA</p>	<p>TZ15.TiO.Si1</p>	

COL SERIES



2-CONNECT®-ABUTMENT SET (STRAIGHT)

incl. 2-CONnect®-Coping and Coping-screw

GH 0,5 mm COL-2CON 8050001S

GH 2,5 mm COL-2CON 8050002S



2-CONNECT®-ABUTMENT

Male Part

GH 0,5 mm COL-2CON 8050001M

GH 2,5 mm COL-2CON 8050002M



2-CONNECT®-SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

KS 61 = 1,8 mm

2CON 9.S3D1.800



2-CONNECT®-COPING

Female Part

2CON F-U



COPING-SCREW

Recommended tightening torques: 25 Ncm

KS 61



2-CONNECT®-DIM ANALOG®

KS 61 = 1,8 mm

2CON 5.DIM.180



SLEEVE (belongs to DIM ANALOG)

plaster model

DIM 5.XSL.A30



2-CONNECT®-IMPLANT PICK-UP

integrated screw

KS 61 = 1,8 mm

2CON TR-00021.8



2-CONNECT®-LAB ANALOG

KS 61 = 1,8 mm

2CON 5.La1.800





2-CONNECT®-ABUTMENT SET (ANGLED)	A: over surface B: over edge	COL-2CON 80500C3S COL-2CON 80500A3S	
2-CONNECT®-ABUTMENT angled 20°, Gingiva Height 2,5 mm Male Part incl. screw, mounting aid	A: over surface B: over edge	COL-2CON 80500C3M COL-2CON 80500A3M	
2-CONNECT®-SCAN BODY 3D-GUIDE® for intra + extra oral scanning KS 62 = 1,6 mm		2CON 9.S3D1.600	
2-CONNECT®-COPING Female Part		2CON F-U	
COPING-SCREW Recommended tightening torques: 15 Ncm		KS 62	
MALE PART SCREW for angled 20° 2-CONnect® Abutment Recommended tightening torques: 20 Ncm		COL 6.Ms5.000	
2-CONNECT®-DIM ANALOG® KS 62 = 1,6 mm		2CON 5.DIM.160	
SLEEVE (belongs to DIM ANALOG) plaster model		DIM 5.XSL.A30	
2-CONNECT®-IMPLANT PICK-UP integrated screw KS 62 = 1,6 mm		2CON TR-00031.6	
2-CONNECT®-LAB ANALOG KS 62 = 1,6 mm		2CON 5.La1.600	