


BEG-SERIES

compatible to:

BEGO Implant Systems®

Semados®

S, RS-RX, RI-Line



PLATFORMS:



3,25 / 3,75 mm



4,1 mm



4,5 mm



5,5 mm



TITANIUM BASE

incl. screw

BEG 80323700



SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

BEG 9.S3D3.237



NT-MILL LAB ANALOG

BEG 5.La3.237



IMPLANT PICK-UP

for open-tray with integrated screw

BEG TR-00013.2



ABUTMENT SCREW



Overview of all screwdrivers on pages 344 and 345

BEG 603255



DIM ANALOG®

BEG 5.DIM.323



SLEEVE (belongs to DIM ANALOG)

plaster model

DIM 5.XSL.A20





<p>NT-PREFORM® incl. coloured milling screw</p>	<p>BEG 10.Mi3.237.05</p>	
<p>SCAN BODY 3D-GUIDE® for intra + extra oral scanning</p>	<p>BEG 9.S3D3.237</p>	
<p>NT-MILL LAB ANALOG for NT-Preform®</p>	<p>BEG 5.La3.237</p>	
<p>COLOURED MILLING SCREW for NT-Preform®</p>	<p>BEG 6.Sc3.255</p>	
<p>ABUTMENT SCREW</p>  <p>Overview of all screwdrivers on pages 344 and 345</p>	<p>BEG 603255</p>	
<p>DIM ANALOG®</p>	<p>BEG 5.DIM.323</p>	
<p>SLEEVE (belongs to DIM ANALOG) plaster model</p>	<p>DIM 5.XSL.A20</p>	
<p>ABUTMENT GRIP for NT-Holder fixture</p>	<p>W 12.100.000</p>	
<p>CUSTOMIZED ABUTMENT milled by NT-iDent</p> <p>Not available in the USA</p>	<p>TZ15.TiO.Si1</p>	



2-CONNECT®-ABUTMENT SET (STRAIGHT)

incl. 2-CONnect®-Coping and Coping-screw

GH 0,5 mm BEG-2CON 8032371S

GH 2,5 mm BEG-2CON 8032372S



2-CONNECT®-ABUTMENT

Male Part

GH 0,5 mm BEG-2CON 8032371M

GH 2,5 mm BEG-2CON 8032372M



2-CONNECT®-SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

KS 61 = 1,8 mm

2CON 9.S3D1.800



2-CONNECT®-COPING

Female Part

2CON F-U



COPING-SCREW

Recommended tightening torques: 25 Ncm

KS 61



2-CONNECT®-DIM ANALOG®

KS 61 = 1,8 mm

2CON 5.DIM.180



SLEEVE (belongs to DIM ANALOG)

plaster model

DIM 5.XSL.A30



2-CONNECT®-IMPLANT PICK-UP

integrated screw

KS 61 = 1,8 mm

2CON TR-00021.8



2-CONNECT®-LAB ANALOG

KS 61 = 1,8 mm

2CON 5.La1.800






2-CONNECT®-ABUTMENT SET (ANGLED)

incl. 2-CONnect®-Coping, mounting aid,
Coping-screw and Male Part Screw, angled 20°,
Gingiva Height 2,5 mm

BEG-2CON 80323A3S



2-CONNECT®-ABUTMENT

angled 20°, Gingiva Height 2,5 mm
Male Part incl. screw,
mounting aid

BEG-2CON 80323A3M



2-CONNECT®-SCAN BODY 3D-GUIDE®

for intra + extra oral scanning
KS 62 = 1,6 mm

2CON 9.S3D1.600



2-CONNECT®-COPING

Female Part

2CON F-U



COPING-SCREW

Recommended tightening torques: 15 Ncm

KS 62



MALE PART SCREW

for angled 20° 2-CONnect® Abutment
Recommended tightening torques: 20 Ncm

BEG 6.Ms3.255



2-CONNECT®-DIM ANALOG®

KS 62 = 1,6 mm

2CON 5.DIM.160



SLEEVE (belongs to DIM ANALOG)

plaster model

DIM 5.XSL.A30



2-CONNECT®-IMPLANT PICK-UP

integrated screw
KS 62 = 1,6 mm

2CON TR-00031.6



2-CONNECT®-LAB ANALOG

KS 62 = 1,6 mm

2CON 5.La1.600





TITANIUM BASE

incl. screw

BEG 80410000



SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

BEG 9.S3D4.100



NT-MILL LAB ANALOG

BEG 5.La4.100



IMPLANT PICK-UP

for open-tray with integrated screw

BEG TR-00024.2



ABUTMENT SCREW



Overview of all screwdrivers on pages 344 and 345

BEG 603255



DIM ANALOG®

BEG 5.DIM.410



SLEEVE (belongs to DIM ANALOG)

plaster model

DIM 5.XSL.A20





<p>NT-PREFORM® incl. coloured milling screw</p>	<p>BEG 10.Mi4.100.05</p>	
<p>SCAN BODY 3D-GUIDE® for intra + extra oral scanning</p>	<p>BEG 9.S3D4.100</p>	
<p>NT-MILL LAB ANALOG for NT-Preform®</p>	<p>BEG 5.La4.100</p>	
<p>COLOURED MILLING SCREW for NT-Preform®</p>	<p>BEG 6.Sc3.255</p>	
<p>ABUTMENT SCREW</p>  <p>Overview of all screwdrivers on pages 344 and 345</p>	<p>BEG 603255</p>	
<p>DIM ANALOG®</p>	<p>BEG 5.DIM.410</p>	
<p>SLEEVE (belongs to DIM ANALOG) plaster model</p>	<p>DIM 5.XSL.A20</p>	
<p>ABUTMENT GRIP for NT-Holder fixture</p>	<p>W 12.200.000</p>	
<p>CUSTOMIZED ABUTMENT milled by NT-iDent</p> <p style="color: orange;">Not available in the USA</p>	<p>TZ15.TiO.Si1</p>	



2-CONNECT®-ABUTMENT SET (STRAIGHT)

incl. 2-CONnect®-Coping and Coping-screw

GH 0,5 mm BEG-2CON 8041001S

GH 2,5 mm BEG-2CON 8041002S



2-CONNECT®-ABUTMENT

Male Part

GH 0,5 mm BEG-2CON 8041001M

GH 2,5 mm BEG-2CON 8041002M



2-CONNECT®-SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

KS 61 = 1,8 mm

2CON 9.S3D1.800



2-CONNECT®-COPING

Female Part

2CON F-U



COPING-SCREW

Recommended tightening torques: 25 Ncm

KS 61



2-CONNECT®-DIM ANALOG®

KS 61 = 1,8 mm

2CON 5.DIM.180



SLEEVE (belongs to DIM ANALOG)

plaster model

DIM 5.XSL.A30



2-CONNECT®-IMPLANT PICK-UP

integrated screw

KS 61 = 1,8 mm

2CON TR-00021.8



2-CONNECT®-LAB ANALOG

KS 61 = 1,8 mm

2CON 5.La1.800






2-CONNECT®-ABUTMENT SET (ANGLED)

incl. 2-CONnect®-Coping, mounting aid,
Coping-screw and Male Part Screw, angled 20°,
Gingiva Height 2,5 mm

BEG-2CON 80410A3S



2-CONNECT®-ABUTMENT

angled 20°, Gingiva Height 2,5 mm
Male Part incl. screw,
mounting aid

BEG-2CON 80410A3M



2-CONNECT®-SCAN BODY 3D-GUIDE®

for intra + extra oral scanning
KS 62 = 1,6 mm

2CON 9.S3D1.600



2-CONNECT®-COPING

Female Part

2CON F-U



COPING-SCREW

Recommended tightening torques: 15 Ncm

KS 62



MALE PART SCREW

for angled 20° 2-CONnect® Abutment
Recommended tightening torques: 20 Ncm

BEG 6.Ms3.255



2-CONNECT®-DIM ANALOG®

KS 62 = 1,6 mm

2CON 5.DIM.160



SLEEVE (belongs to DIM ANALOG)

plaster model

DIM 5.XSL.A30



2-CONNECT®-IMPLANT PICK-UP

integrated screw
KS 62 = 1,6 mm

2CON TR-00031.6



2-CONNECT®-LAB ANALOG

KS 62 = 1,6 mm

2CON 5.La1.600





TITANIUM BASE

incl. screw

BEG 80450000



SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

BEG 9.S3D4.500



NT-MILL LAB ANALOG

BEG 5.La4.500



IMPLANT PICK-UP

for open-tray with integrated screw

BEG TR-00034.5



ABUTMENT SCREW



Overview of all screwdrivers on pages 344 and 345

BEG 603255



DIM ANALOG®

BEG 5.DIM.450



SLEEVE (belongs to DIM ANALOG)

plaster model

DIM 5.XSL.A30





<p>NT-PREFORM® incl. coloured milling screw</p>	<p>BEG 10.Mi4.500.05</p>	
<p>SCAN BODY 3D-GUIDE® for intra + extra oral scanning</p>	<p>BEG 9.S3D4.500</p>	
<p>NT-MILL LAB ANALOG for NT-Preform®</p>	<p>BEG 5.La4.500</p>	
<p>COLOURED MILLING SCREW for NT-Preform®</p>	<p>BEG 6.Sc3.255</p>	
<p>ABUTMENT SCREW</p>  <p>Overview of all screwdrivers on pages 344 and 345</p>	<p>BEG 603255</p>	
<p>DIM ANALOG®</p>	<p>BEG 5.DIM.450</p>	
<p>SLEEVE (belongs to DIM ANALOG) plaster model</p>	<p>DIM 5.XSL.A30</p>	
<p>ABUTMENT GRIP for NT-Holder fixture</p>	<p>W 12.300.000</p>	
<p>CUSTOMIZED ABUTMENT milled by NT-iDent</p> <p>Not available in the USA</p>	<p>TZ15.TiO.Si1</p>	



2-CONNECT®-ABUTMENT SET (STRAIGHT)

incl. 2-CONnect®-Coping and Coping-screw

GH 0,5 mm BEG-2CON 8045001S

GH 2,5 mm BEG-2CON 8045002S



2-CONNECT®-ABUTMENT

Male Part

GH 0,5 mm BEG-2CON 8045001M

GH 2,5 mm BEG-2CON 8045002M



2-CONNECT®-SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

KS 61 = 1,8 mm

2CON 9.S3D1.800



2-CONNECT®-COPING

Female Part

2CON F-U



COPING-SCREW

Recommended tightening torques: 25 Ncm

KS 61



2-CONNECT®-DIM ANALOG®

KS 61 = 1,8 mm

2CON 5.DIM.180



SLEEVE (belongs to DIM ANALOG)

plaster model

DIM 5.XSL.A30



2-CONNECT®-IMPLANT PICK-UP

integrated screw

KS 61 = 1,8 mm

2CON TR-00021.8



2-CONNECT®-LAB ANALOG

KS 61 = 1,8 mm

2CON 5.La1.800






2-CONNECT®-ABUTMENT SET (ANGLED)

incl. 2-CONnect®-Coping, mounting aid,
Coping-screw and Male Part Screw, angled 20°,
Gingiva Height 2,5 mm

BEG-2CON 80450A3S



2-CONNECT®-ABUTMENT

angled 20°, Gingiva Height 2,5 mm
Male Part incl. screw,
mounting aid

BEG-2CON 80450A3M



2-CONNECT®-SCAN BODY 3D-GUIDE®

for intra + extra oral scanning
KS 62 = 1,6 mm

2CON 9.S3D1.600



2-CONNECT®-COPING

Female Part

2CON F-U



COPING-SCREW

Recommended tightening torques: 15 Ncm

KS 62



MALE PART SCREW

for angled 20° 2-CONnect® Abutment
Recommended tightening torques: 20 Ncm

BEG 6.Ms3.255



2-CONNECT®-DIM ANALOG®

KS 62 = 1,6 mm

2CON 5.DIM.160



SLEEVE (belongs to DIM ANALOG)

plaster model

DIM 5.XSL.A30



2-CONNECT®-IMPLANT PICK-UP

integrated screw
KS 62 = 1,6 mm

2CON TR-00031.6



2-CONNECT®-LAB ANALOG

KS 62 = 1,6 mm

2CON 5.La1.600





TITANIUM BASE

incl. screw

BEG 80550000



SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

BEG 9.S3D5.500



NT-MILL LAB ANALOG

BEG 5.La5.500



IMPLANT PICK-UP

for open-tray with integrated screw

BEG TR-00045.5



ABUTMENT SCREW



Overview of all screwdrivers on pages 344 and 345

BEG 603255



DIM ANALOG®

BEG 5.DIM.550



SLEEVE (belongs to DIM ANALOG)


plaster model

DIM 5.XSL.A30





<p>NT-PREFORM® incl. coloured milling screw</p>	<p>BEG 10.Mi5.500.05</p>	
<p>SCAN BODY 3D-GUIDE® for intra + extra oral scanning</p>	<p>BEG 9.S3D5.500</p>	
<p>NT-MILL LAB ANALOG for NT-Preform®</p>	<p>BEG 5.La5.500</p>	 groove
<p>COLOURED MILLING SCREW for NT-Preform®</p>	<p>BEG 6.Sc3.255</p>	
<p>ABUTMENT SCREW</p>  <p>Overview of all screwdrivers on pages 344 and 345</p>	<p>BEG 603255</p>	
<p>DIM ANALOG®</p>	<p>BEG 5.DIM.550</p>	
<p>SLEEVE (belongs to DIM ANALOG) plaster model</p>	<p>DIM 5.XSL.A30</p>	
<p>ABUTMENT GRIP for NT-Holder fixture</p>	<p>W 12.300.000</p>	
<p>CUSTOMIZED ABUTMENT milled by NT-iDent</p> <p style="color: red;">Not available in the USA</p>	<p>TZ15.TiO.Si1</p>	



2-CONNECT®-ABUTMENT SET (STRAIGHT)

incl. 2-CONnect®-Coping and Coping-screw

GH 0,5 mm BEG-2CON 8055001S

GH 2,5 mm BEG-2CON 8055002S



2-CONNECT®-ABUTMENT

Male Part

GH 0,5 mm BEG-2CON 8055001M

GH 2,5 mm BEG-2CON 8055002M



2-CONNECT®-SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

KS 61 = 1,8 mm

2CON 9.S3D1.800



2-CONNECT®-COPING

Female Part

2CON F-U



COPING-SCREW

Recommended tightening torques: 25 Ncm

KS 61



2-CONNECT®-DIM ANALOG®

KS 61 = 1,8 mm

2CON 5.DIM.180



SLEEVE (belongs to DIM ANALOG)

plaster model

DIM 5.XSL.A30



2-CONNECT®-IMPLANT PICK-UP

integrated screw

KS 61 = 1,8 mm

2CON TR-00021.8



2-CONNECT®-LAB ANALOG

KS 61 = 1,8 mm

2CON 5.La1.800






2-CONNECT®-ABUTMENT SET (ANGLED)

incl. 2-CONnect®-Coping, mounting aid,
Coping-screw and Male Part Screw, angled 20°,
Gingiva Height 2,5 mm

BEG-2CON 80550A3S



2-CONNECT®-ABUTMENT

angled 20°, Gingiva Height 2,5 mm
Male Part incl. screw,
mounting aid

BEG-2CON 80550A3M



2-CONNECT®-SCAN BODY 3D-GUIDE®

for intra + extra oral scanning
KS 62 = 1,6 mm

2CON 9.S3D1.600



2-CONNECT®-COPING

Female Part

2CON F-U



COPING-SCREW

Recommended tightening torques: 15 Ncm

KS 62



MALE PART SCREW

for angled 20° 2-CONnect® Abutment
Recommended tightening torques: 20 Ncm

BEG 6.Ms3.255



2-CONNECT®-DIM ANALOG®

KS 62 = 1,6 mm

2CON 5.DIM.160



SLEEVE (belongs to DIM ANALOG)

plaster model

DIM 5.XSL.A30



2-CONNECT®-IMPLANT PICK-UP

integrated screw
KS 62 = 1,6 mm

2CON TR-00031.6



2-CONNECT®-LAB ANALOG

KS 62 = 1,6 mm

2CON 5.La1.600

