

H-SERIES

compatible to:

Biomet 3i[®]

Certain[®]



PLATFORMS:



3,4 mm NP



4,1 mm RP



5,0 mm WP



TITANIUM BASE

incl. screw

H 800



SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

H 9.S3D3.400



LAB ANALOG

H 50



IMPLANT PICK-UP

for open-tray with integrated screw

H TR-NP013.4



ABUTMENT SCREW



Overview of all screwdrivers on pages 344 and 345

H 60



DIM ANALOG®

H 5.DIM.340



SLEEVE (belongs to DIM ANALOG)

plaster model

DIM 5.XSL.A10



H
SERIES



<p>NT-PREFORM® incl. coloured milling screw</p>	<p>GH 0,5 mm GH 1,3 mm</p>	<p>H10.Mi3.400.05 H10.Mi3.400</p>	
<p>SCAN BODY 3D-GUIDE® for intra + extra oral scanning</p>		<p>H 9.S3D3.400</p>	
<p>NT-MILL LAB ANALOG for NT-Preform®</p>		<p>H 5.La3.400</p>	
<p>COLOURED MILLING SCREW for NT-Preform®</p>		<p>H 6.Sc3.450</p>	
<p>ABUTMENT SCREW</p>	<p>Overview of all screwdrivers on pages 344 and 345</p>	<p>H 60</p>	
<p>DIM ANALOG®</p>		<p>H 5.DIM.340</p>	
<p>SLEEVE (belongs to DIM ANALOG®) plaster model</p>		<p>DIM 5.XSL.A10</p>	
<p>ABUTMENT GRIP for NT-Holder fixture</p>		<p>W 12.100.000</p>	
<p>CUSTOMIZED ABUTMENT milled by NT-iDent</p>	<p>Not available in the USA</p>	<p>TZ15.TiO.Si1</p>	

H
SERIES



2-CONNECT®-ABUTMENT SET (STRAIGHT)

incl. 2-CONnect®-Coping and Coping-screw

GH 0,5 mm

H 800 S

GH 2,5 mm

H-2CON 8034002S



2-CONNECT®-ABUTMENT

Male Part

GH 0,5 mm

H 800 M

GH 2,5 mm

H-2CON 8034002M



2-CONNECT®-SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

KS 60 = 1,4 mm

2CON 9.S3D1.400



2-CONNECT®-COPING

Female Part

2CON F-U



COPING-SCREW

Recommended tightening torques: 15 Ncm

KS 60



2-CONNECT®-DIM ANALOG®

KS 60 = 1,4 mm

2CON 5.DIM.140



SLEEVE (belongs to DIM ANALOG®)

plaster model

DIM 5.XSL.A30



2-CONNECT®-IMPLANT PICK-UP

integrated screw

KS 60 = 1,4 mm

CON TR-00011.4



2-CONNECT®-LAB ANALOG

KS 60 = 1,4 mm

2CON 5.La1.400





<p>2-CONNECT®-ABUTMENT SET (ANGLED) incl. 2-CONnect®-Coping, mounting aid, Coping-screw and Male Part Screw, angled 20°, Gingiva Height 2,5 mm</p>	<p>H-2CON 80340A3S</p>	
<p>2-CONNECT®-ABUTMENT angled 20°, Gingiva Height 2,5 mm Male Part incl. screw, mounting aid</p>	<p>H-2CON 80340A3M</p>	
<p>2-CONNECT®-SCAN BODY 3D-GUIDE® for intra + extra oral scanning KS 62 = 1,6 mm</p>	<p>2CON 9.S3D1.600</p>	
<p>2-CONNECT®-COPING Female Part</p>	<p>2CON F-U</p>	
<p>COPING-SCREW Recommended tightening torques: 15 Ncm</p>	<p>KS 62</p>	
<p>MALE PART SCREW for angled 20° 2-CONnect® Abutment Recommended tightening torques: 20 Ncm</p>	<p>H 6.Ms3.450</p>	
<p>2-CONNECT®-DIM ANALOG® KS 62 = 1,6 mm</p>	<p>2CON 5.DIM.160</p>	
<p>SLEEVE (belongs to DIM ANALOG®) plaster model</p>	<p>DIM 5.XSL.A30</p>	
<p>2-CONNECT®-IMPLANT PICK-UP integrated screw KS 62 = 1,6 mm</p>	<p>2CON TR-00031.6</p>	
<p>2-CONNECT®-LAB ANALOG KS 62 = 1,6 mm</p>	<p>2CON 5.La1.600</p>	

H
SERIES



TITANIUM BASE

incl. screw

H 810



SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

H 9.S3D4.150



LAB ANALOG

H 51



IMPLANT PICK-UP

for open-tray with integrated screw

H TR-RP024.1



ABUTMENT SCREW



Overview of all screwdrivers on pages 344 and 345

H 60



DIM ANALOG®

H 5.DIM.410



SLEEVE (belongs to DIM ANALOG)

plaster model

DIM 5.XSL.A20



<p>NT-PREFORM® incl. coloured milling screw</p>	<p>GH 0,5 mm GH 1,3 mm</p>	<p>H 10.Mi4.100.05 H 10.Mi4.100</p>	
<p>SCAN BODY 3D-GUIDE® for intra + extra oral scanning</p>		<p>H 9.S3D4.150</p>	
<p>NT-MILL LAB ANALOG for NT-Preform®</p>		<p>H 5.La4.100</p>	 groove
<p>COLOURED MILLING SCREW for NT-Preform®</p>		<p>H 6.Sc3.450</p>	
<p>ABUTMENT SCREW</p> 	<p>Overview of all screwdrivers on pages 344 and 345</p>	<p>H 60</p>	
<p>DIM ANALOG®</p>		<p>H 5.DIM.410</p>	
<p>SLEEVE (belongs to DIM ANALOG®) plaster model</p>		<p>DIM 5.XSL.A20</p>	
<p>ABUTMENT GRIP for NT-Holder fixture</p>		<p>W 12.200.000</p>	
<p>CUSTOMIZED ABUTMENT milled by NT-iDent</p>	<p>Not available in the USA</p>	<p>TZ15.TiO.Si1</p>	

H
SERIES



2-CONNECT®-ABUTMENT SET (STRAIGHT)

incl. 2-CONnect®-Coping and Coping-screw

GH 0,5 mm

H 810 S

GH 2,5 mm

H-2CON 8041002S



2-CONNECT®-ABUTMENT

Male Part

GH 0,5 mm

H 810 M

GH 2,5 mm

H-2CON 8041002M



2-CONNECT®-SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

KS 60 = 1,4 mm

2CON 9.S3D1.400



2-CONNECT®-COPING

Female Part

2CON F-U



COPING-SCREW

Recommended tightening torques: 15 Ncm

KS 60



2-CONNECT®-DIM ANALOG®

KS 60 = 1,4 mm

2CON 5.DIM.140



SLEEVE (belongs to DIM ANALOG®)

plaster model

DIM 5.XSL.A30



2-CONNECT®-IMPLANT PICK-UP

integrated screw

KS 60 = 1,4 mm

2CON TR-00011.4



2-CONNECT®-LAB ANALOG

KS 60 = 1,4 mm

2CON 5.La1.400





<p>2-CONNECT®-ABUTMENT SET (ANGLED) incl. 2-CONnect®-Coping, mounting aid, Coping-screw and Male Part Screw, angled 20°, Gingiva Height 2,5 mm</p>	<p>H-2CON 80410A3S</p>	
<p>2-CONNECT®-ABUTMENT angled 20°, Gingiva Height 2,5 mm Male Part incl. screw, mounting aid</p>	<p>H-2CON 80410A3M</p>	
<p>2-CONNECT®-SCAN BODY 3D-GUIDE® for intra + extra oral scanning KS 62 = 1,6 mm</p>	<p>2CON 9.S3D1.600</p>	
<p>2-CONNECT®-COPING Female Part</p>	<p>2CON F-U</p>	
<p>COPING-SCREW Recommended tightening torques: 15 Ncm</p>	<p>KS 62</p>	
<p>MALE PART SCREW for angled 20° 2-CONnect® Abutment Recommended tightening torques: 20 Ncm</p>	<p>H 6.Ms3.450</p>	
<p>2-CONNECT®-DIM ANALOG® KS 62 = 1,6 mm</p>	<p>2CON 5.DIM.160</p>	
<p>SLEEVE (belongs to DIM ANALOG®) plaster model</p>	<p>DIM 5.XSL.A30</p>	
<p>2-CONNECT®-IMPLANT PICK-UP integrated screw KS 62 = 1,6 mm</p>	<p>2CON TR-00031.6</p>	
<p>2-CONNECT®-LAB ANALOG KS 62 = 1,6 mm</p>	<p>2CON 5.La1.600</p>	



TITANIUM BASE

incl. screw

H 820



SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

H 9.S3D4.150



LAB ANALOG

H 52



IMPLANT PICK-UP

for open-tray with integrated screw

H TR-WP035.0



ABUTMENT SCREW



Overview of all screwdrivers on pages 344 and 345

H 60



DIM ANALOG®

H 5.DIM.500



SLEEVE (belongs to DIM ANALOG)

plaster model

DIM 5.XSL.A30





<p>NT-PREFORM® incl. coloured milling screw</p>	<table border="1"> <tr> <td data-bbox="798 403 933 448">GH 0,5 mm</td> <td data-bbox="933 403 1220 448">H 10.Mi5.000.05</td> </tr> <tr> <td data-bbox="798 448 933 492">GH 1,3 mm</td> <td data-bbox="933 448 1220 492">H 10.Mi5.000</td> </tr> </table>	GH 0,5 mm	H 10.Mi5.000.05	GH 1,3 mm	H 10.Mi5.000	
GH 0,5 mm	H 10.Mi5.000.05					
GH 1,3 mm	H 10.Mi5.000					
<p>SCAN BODY 3D-GUIDE® for intra + extra oral scanning</p>	<p>H 9.S3D4.150</p>					
<p>NT-MILL LAB ANALOG for NT-Preform®</p>	<p>H 5.La5.000</p>					
<p>COLOURED MILLING SCREW for NT-Preform®</p>	<p>H 6.Sc3.450</p>					
<p>ABUTMENT SCREW</p>	<p>H 60</p>					
<p>Overview of all screwdrivers on pages 344 and 345</p>						
<p>DIM ANALOG®</p>	<p>H 5.DIM.500</p>					
<p>SLEEVE (belongs to DIM ANALOG®) plaster model</p>	<p>DIM 5.XSL.A30</p>					
<p>ABUTMENT GRIP for NT-Holder fixture</p>	<p>W 12.300.000</p>					
<p>CUSTOMIZED ABUTMENT milled by NT-iDent</p>	<p>Not available in the USA TZ15.TiO.Si1</p>					

H
SERIES