

T-SERIES

compatible to:

DENTSPLY Implants[®]
Xive[®]



PLATFORMS:



3,4 mm



3,8 mm



4,5 mm



5,5 mm



TITANIUM BASE

incl. screw

T 800



SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

T 9.S3D3.400



LAB ANALOG

T 50



IMPLANT PICK-UP

for open-tray with integrated screw

T TR-00013.4



ABUTMENT SCREW



Overview of all screwdrivers on pages 344 and 345

T 60



DIM ANALOG®

T 5.DIM.340



SLEEVE (belongs to DIM ANALOG)

plaster model

DIM 5.XSL.A10





<p>NT-PREFORM® incl. coloured milling screw</p>	<p>GH 0,5 mm T 10.Mi3.400.05 GH 1,3 mm T 10.Mi3.400</p>	
<p>SCAN BODY 3D-GUIDE® for intra + extra oral scanning</p>	<p>T 9.S3D3.400</p>	
<p>NT-MILL LAB ANALOG for NT-Preform®</p>	<p>T 5.La3.400</p>	
<p>COLOURED MILLING SCREW for NT-Preform®</p>	<p>T 6.Sc3.455</p>	
<p>ABUTMENT SCREW</p>  <p>Overview of all screwdrivers on pages 344 and 345</p>	<p>T 60</p>	
<p>DIM ANALOG®</p>	<p>T 5.DIM.340</p>	
<p>SLEEVE (belongs to DIM ANALOG®) plaster model</p>	<p>DIM 5.XSL.A10</p>	
<p>ABUTMENT GRIP for NT-Holder fixture</p>	<p>W 12.100.000</p>	
<p>CUSTOMIZED ABUTMENT milled by NT-iDent</p>	<p>Not available in the USA T Z15.TiO.Si1</p>	



2-CONNECT®-ABUTMENT SET (STRAIGHT)

incl. 2-CONnect®-Coping and Coping-screw

GH 0,5 mm	T 800 S
GH 2,5 mm	T-2CON 8034002S



2-CONNECT®-ABUTMENT

Male Part

GH 0,5 mm	T 800 M
GH 2,5 mm	T-2CON 8034002M



2-CONNECT®-SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

KS 60 = 1,4 mm

2CON 9.S3D1.400



2-CONNECT®-COPING

Female Part

2CON F-U



COPING-SCREW

Recommended tightening torques: 15 Ncm

KS 60



2-CONNECT®-DIM ANALOG®

KS 60 = 1,4 mm

2CON 5.DIM.140



SLEEVE (belongs to DIM ANALOG®)

plaster model

DIM 5.XSL.A30



2-CONNECT®-IMPLANT PICK-UP

integrated screw

KS 60 = 1,4 mm

2CON TR-00011.4



2-CONNECT®-LAB ANALOG

KS 60 = 1,4 mm

2CON 5.La1.400





<p>2-CONNECT®-ABUTMENT SET (ANGLED) incl. 2-CONnect®-Coping, mounting aid, Coping-screw and Male Part Screw, angled 20°, Gingiva Height 2,5 mm</p>	<p>T-2CON 80340A3S</p>	
<p>2-CONNECT®-ABUTMENT angled 20°, Gingiva Height 2,5 mm Male Part incl. screw, mounting aid</p>	<p>T-2CON 80340A3M</p>	
<p>2-CONNECT®-SCAN BODY 3D-GUIDE® for intra + extra oral scanning KS 62 = 1,6 mm</p>	<p>2CON 9.S3D1.600</p>	
<p>2-CONNECT®-COPING Female Part</p>	<p>2CON F-U</p>	
<p>COPING-SCREW Recommended tightening torques: 15 Ncm</p>	<p>KS 62</p>	
<p>MALE PART SCREW for angled 20° 2-CONnect® Abutment Recommended tightening torques: 20 Ncm</p>	<p>T 6.Ms3.455</p>	
<p>2-CONNECT®-DIM ANALOG® KS 62 = 1,6 mm</p>	<p>2CON 5.DIM.160</p>	
<p>SLEEVE (belongs to DIM ANALOG®) plaster model</p>	<p>DIM 5.XSL.A30</p>	
<p>2-CONNECT®-IMPLANT PICK-UP integrated screw KS 62 = 1,6 mm</p>	<p>2CON TR-00031.6</p>	
<p>2-CONNECT®-LAB ANALOG KS 62 = 1,6 mm</p>	<p>2CON 5.La1.600</p>	



TITANIUM BASE

incl. screw

T 805



SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

T 9.S3D3.800



LAB ANALOG

T 55



IMPLANT PICK-UP

for open-tray with integrated screw

T TR-00023.8



ABUTMENT SCREW



Overview of all screwdrivers on pages 344 and 345

T 60



DIM ANALOG®

T 5.DIM.380



SLEEVE (belongs to DIM ANALOG)

plaster model

DIM 5.XSL.A20



<p>NT-PREFORM® incl. coloured milling screw</p>	<p>GH 0,5 mm GH 1,3 mm</p>	<p>T 10.Mi3.800.05 T 10.Mi3.800</p>	
<p>SCAN BODY 3D-GUIDE® for intra + extra oral scanning</p>		<p>T 9.S3D3.800</p>	
<p>NT-MILL LAB ANALOG for NT-Preform®</p>		<p>T 5.La3.800</p>	 groove
<p>COLOURED MILLING SCREW for NT-Preform®</p>		<p>T 6.Sc3.455</p>	
<p>ABUTMENT SCREW</p>  <p>Overview of all screwdrivers on pages 344 and 345</p>		<p>T 60</p>	
<p>DIM ANALOG®</p>		<p>T 5.DIM.380</p>	
<p>SLEEVE (belongs to DIM ANALOG®) plaster model</p>		<p>DIM 5.XSL.A20</p>	
<p>ABUTMENT GRIP for NT-Holder fixture</p>		<p>W 12.200.000</p>	
<p>CUSTOMIZED ABUTMENT milled by NT-iDent</p>	<p>Not available in the USA</p>	<p>TZ15.TiO.Si1</p>	



2-CONNECT®-ABUTMENT SET (STRAIGHT)

incl. 2-CONnect®-Coping and Coping-screw

GH 0,5 mm	T 805 S
GH 2,5 mm	T-2CON 8038002S



2-CONNECT®-ABUTMENT

Male Part

GH 0,5 mm	T 805 M
GH 2,5 mm	T-2CON 8038002M



2-CONNECT®-SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

KS 61 = 1,8 mm

2CON 9.S3D1.800



2-CONNECT®-COPING

Female Part

2CON F-U



COPING-SCREW

Recommended tightening torques: 20 Ncm

KS 61



2-CONNECT®-DIM ANALOG®

KS 61 = 1,8 mm

2CON 5.DIM.180



SLEEVE (belongs to DIM ANALOG®)

plaster model

DIM 5.XSL.A30



2-CONNECT®-IMPLANT PICK-UP

integrated screw

KS 61 = 1,8 mm

2CON TR-00021.8



2-CONNECT®-LAB ANALOG

KS 61 = 1,8 mm

2CON 5.La1.800





<p>2-CONNECT®-ABUTMENT SET (ANGLED) incl. 2-CONnect®-Coping, mounting aid, Coping-screw and Male Part Screw, angled 20°, Gingiva Height 2,5 mm</p>	T-2CON 80380A3S	
<p>2-CONNECT®-ABUTMENT angled 20°, Gingiva Height 2,5 mm Male Part incl. screw, mounting aid</p>	T-2CON 80380A3M	
<p>2-CONNECT®-SCAN BODY 3D-GUIDE® for intra + extra oral scanning KS 62 = 1,6 mm</p>	2CON 9.S3D1.600	
<p>2-CONNECT®-COPING Female Part</p>	2CON F-U	
<p>COPING-SCREW Recommended tightening torques: 15 Ncm</p>	KS 62	
<p>MALE PART SCREW for angled 20° 2-CONnect® Abutment Recommended tightening torques: 20 Ncm</p>	T 6.Ms3.455	
<p>2-CONNECT®-DIM ANALOG® KS 62 = 1,6 mm</p>	2CON 5.DIM.160	
<p>SLEEVE (belongs to DIM ANALOG®) plaster model</p>	DIM 5.XSL.A30	
<p>2-CONNECT®-IMPLANT PICK-UP integrated screw KS 62 = 1,6 mm</p>	2CON TR-00031.6	
<p>2-CONNECT®-LAB ANALOG KS 62 = 1,6 mm</p>	2CON 5.La1.600	



TITANIUM BASE

incl. screw

T 810



SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

T 9.S3D4.555



LAB ANALOG

T 51



IMPLANT PICK-UP

for open-tray with integrated screw

T TR-00034.5



ABUTMENT SCREW



Overview of all screwdrivers on pages 344 and 345

T 60



DIM ANALOG®

T 5.DIM.450



SLEEVE (belongs to DIM ANALOG)

plaster model

DIM 5.XSL.A20





NT-PREFORM®	incl. coloured milling screw	GH 0,5 mm	T 10.Mi4.500.05	
		GH 1,3 mm	T 10.Mi4.500	
SCAN BODY 3D-GUIDE®	for intra + extra oral scanning		T 9.S3D4.555	
NT-MILL LAB ANALOG	for NT-Preform®		T 5.La4.500	
COLOURED MILLING SCREW	for NT-Preform®		T 6.Sc3.455	
ABUTMENT SCREW		Overview of all screwdrivers on pages 344 and 345	T 60	
DIM ANALOG®			T 5.DIM.450	
SLEEVE (belongs to DIM ANALOG®)	plaster model		DIM 5.XSL.A20	
ABUTMENT GRIP	for NT-Holder fixture		W 12.200.000	
CUSTOMIZED ABUTMENT	milled by NT-iDent	Not available in the USA	TZ15.TiO.Si1	



2-CONNECT®-ABUTMENT SET (STRAIGHT)

incl. 2-CONnect®-Coping and Coping-screw

GH 0,5 mm

T 810 S

GH 2,5 mm

T-2CON 8045002S



2-CONNECT®-ABUTMENT

Male Part

GH 0,5 mm

T 810 M

GH 2,5 mm

T-2CON 8045002M



2-CONNECT®-SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

KS 61 = 1,8 mm

2CON 9.S3D1.800



2-CONNECT®-COPING

Female Part

2CON F-U



COPING-SCREW

Recommended tightening torques: 20 Ncm

KS 61



2-CONNECT®-DIM ANALOG®

KS 61 = 1,8 mm

2CON 5.DIM.180



SLEEVE (belongs to DIM ANALOG®)

plaster model

DIM 5.XSL.A30



2-CONNECT®-IMPLANT PICK-UP

integrated screw

KS 61 = 1,8 mm

2CON TR-00021.8



2-CONNECT®-LAB ANALOG

KS 61 = 1,8 mm

2CON 5.La1.800





<p>2-CONNECT®-ABUTMENT SET (ANGLED) incl. 2-CONnect®-Coping, mounting aid, Coping-screw and Male Part Screw, angled 20°, Gingiva Height 2,5 mm</p>	<p>T-2CON 80450A3S</p>	
<p>2-CONNECT®-ABUTMENT angled 20°, Gingiva Height 2,5 mm Male Part incl. screw, mounting aid</p>	<p>T-2CON 80450A3M</p>	
<p>2-CONNECT®-SCAN BODY 3D-GUIDE® for intra + extra oral scanning KS 62 = 1,6 mm</p>	<p>2CON 9.S3D1.600</p>	
<p>2-CONNECT®-COPING Female Part</p>	<p>2CON F-U</p>	
<p>COPING-SCREW Recommended tightening torques: 15 Ncm</p>	<p>KS 62</p>	
<p>MALE PART SCREW for angled 20° 2-CONnect® Abutment Recommended tightening torques: 20 Ncm</p>	<p>T 6.Ms3.455</p>	
<p>2-CONNECT®-DIM ANALOG® KS 62 = 1,6 mm</p>	<p>2CON 5.DIM.160</p>	
<p>SLEEVE (belongs to DIM ANALOG®) plaster model</p>	<p>DIM 5.XSL.A30</p>	
<p>2-CONNECT®-IMPLANT PICK-UP integrated screw KS 62 = 1,6 mm</p>	<p>2CON TR-00031.6</p>	
<p>2-CONNECT®-LAB ANALOG KS 62 = 1,6 mm</p>	<p>2CON 5.La1.600</p>	



TITANIUM BASE

incl. screw

T 820



SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

T 9.S3D4.555



LAB ANALOG

T 52



IMPLANT PICK-UP

for open-tray with integrated screw

T TR-00034.5



ABUTMENT SCREW



Overview of all screwdrivers on pages 344 and 345

T 60



DIM ANALOG®

T 5.DIM.550



SLEEVE (belongs to DIM ANALOG)

plaster model

DIM 5.XSL.A30





NT-PREFORM®	incl. coloured milling screw	GH 0,5 mm	T 10.Mi5.500.05	
		GH 1,3 mm	T 10.Mi5.500	
SCAN BODY 3D-GUIDE®	for intra + extra oral scanning		T 9.S3D4.555	
NT-MILL LAB ANALOG	for NT-Preform®		T 5.La5.500	
COLOURED MILLING SCREW	for NT-Preform®		T 6.Sc3.455	
ABUTMENT SCREW		Overview of all screwdrivers on pages 344 and 345	T 60	
DIM ANALOG®			T 5.DIM.550	
SLEEVE (belongs to DIM ANALOG®)	plaster model		DIM 5.XSL.A30	
ABUTMENT GRIP	for NT-Holder fixture		W 12.300.000	
CUSTOMIZED ABUTMENT	milled by NT-iDent	Not available in the USA	TZ15.TiO.Si1	