


# MEX-SERIES

compatible to:

*Medentis Medical*<sup>®</sup>  
*ICX*<sup>®</sup>



**PLATFORMS:**



*3,7 - 4,8 mm*



## TITANIUM BASE

incl. screw

MEX 80374800



## SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

MEX 9.S3D3.748



## NT-MILL LAB ANALOG

MEX 5.La3.748



## IMPLANT PICK-UP

for open-tray with integrated screw

MEX TR-00013.7



## ABUTMENT SCREW



Overview of all screwdrivers on pages 344 and 345

MEX 603748



## DIM ANALOG®

MEX 5.DIM.374



## SLEEVE (belongs to DIM ANALOG)

plaster model

DIM 5.XSL.A20





<p><b>NT-PREFORM®</b> incl. coloured milling screw</p>	<p>MEX 10.Mi3.748.05</p>	
<p><b>SCAN BODY 3D-GUIDE®</b> for intra + extra oral scanning</p>	<p>MEX 9.S3D3.748</p>	
<p><b>NT-MILL LAB ANALOG</b> for NT-Preform®</p>	<p>MEX 5.La3.748</p>	
<p><b>COLOURED MILLING SCREW</b> for NT-Preform®</p>	<p>MEX 6.Sc3.748</p>	
<p><b>ABUTMENT SCREW</b></p>  <p>Overview of all screwdrivers on pages 344 and 345</p>	<p>MEX 603748</p>	
<p><b>DIM ANALOG®</b></p>	<p>MEX 5.DIM.374</p>	
<p><b>SLEEVE</b> (belongs to DIM ANALOG®) plaster model</p>	<p>DIM 5.XSL.A20</p>	
<p><b>ABUTMENT GRIP</b> for NT-Holder fixture</p>	<p>W 12.200.000</p>	
<p><b>CUSTOMIZED ABUTMENT</b> milled by NT-iDent</p> <p>Not available in the USA</p>	<p>TZ15.TiO.Si1</p>	

**MEX**  
SERIES



### 2-CONNECT®-ABUTMENT SET (STRAIGHT)

incl. 2-CONnect®-Coping and Coping-screw

GH 0,5 mm MEX-2CON 8037481S

GH 2,5 mm MEX-2CON 8037482S



### 2-CONNECT®-ABUTMENT

Male Part

GH 0,5 mm MEX-2CON 8037481M

GH 2,5 mm MEX-2CON 8037482M



### 2-CONNECT®-SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

KS 61 = 1,8 mm

2CON 9.S3D1.800



### 2-CONNECT®-COPING

Female Part

2CON F-U



### COPING-SCREW

Recommended tightening torques: 25 Ncm

KS 61



### 2-CONNECT®-DIM ANALOG®

KS 61 = 1,8 mm

2CON 5.DIM.180



### SLEEVE (belongs to DIM ANALOG®)

plaster model

DIM 5.XSL.A30



### 2-CONNECT®-IMPLANT PICK-UP

integrated screw

KS 61 = 1,8 mm

2CON TR-00021.8



### 2-CONNECT®-LAB ANALOG

KS 61 = 1,8 mm

2CON 5.La1.800





<p><b>2-CONNECT®-ABUTMENT SET (ANGLED)</b> incl. 2-CONnect®-Coping, mounting aid, Coping-screw and Male Part Screw, angled 20°, Gingiva Height 2,5 mm</p>	<p>MEX-2CON 80374A3S</p>	
<p><b>2-CONNECT®-ABUTMENT</b> angled 20°, Gingiva Height 2,5 mm Male Part incl. screw, mounting aid</p>	<p>MEX-2CON 80374A3M</p>	
<p><b>2-CONNECT®-SCAN BODY 3D-GUIDE®</b> for intra + extra oral scanning KS 62 = 1,6 mm</p>	<p>2CON 9.S3D1.600</p>	
<p><b>2-CONNECT®-COPING</b> Female Part</p>	<p>2CON F-U</p>	
<p><b>COPING-SCREW</b> Recommended tightening torques: 15 Ncm</p>	<p>KS 62</p>	
<p><b>MALE PART SCREW</b> for angled 20° 2-CONnect® Abutment Recommended tightening torques: 20 Ncm</p>	<p>MEX 6.Ms3.748</p>	
<p><b>2-CONNECT®-DIM ANALOG®</b> KS 62 = 1,6 mm</p>	<p>2CON 5.DIM.160</p>	
<p><b>SLEEVE</b> (belongs to DIM ANALOG®) plaster model</p>	<p>DIM 5.XSL.A30</p>	
<p><b>2-CONNECT®-IMPLANT PICK-UP</b> integrated screw KS 62 = 1,6 mm</p>	<p>2CON TR-00031.6</p>	
<p><b>2-CONNECT®-LAB ANALOG</b> KS 62 = 1,6 mm</p>	<p>2CON 5.La1.600</p>	