


MEG-SERIES

compatible to:

Megagen[®]
AnyRidge[®]



PLATFORM:



3,5 - 5,5 mm



TITANIUM BASE

incl. screw

MEG 80355500



SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

MEG 9.S3D3.555



NT-MILL LAB ANALOG

MEG 5.La3.555



IMPLANT PICK-UP

for open-tray with integrated screw

MEG TR-00013.5



ABUTMENT SCREW



Overview of all screwdrivers on pages 344 and 345 in the catalogue

MEG 603555



DIM ANALOG®

MEG 5.DIM.355



SLEEVE (belongs to DIM ANALOG)

plaster model

DIM 5.XSL.A20



MEG-SERIES

TITANIUM BASE NX

incl. screw

GH 0,3mm

MEG 82355003

GH 1,3mm

MEG 82355013



SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

MEG 9.S3D3.555



NT-MILL LAB ANALOG

MEG 5.La3.555



IMPLANT PICK-UP

for open-tray with integrated screw

MEG TR-00013.5



ABUTMENT SCREW



Overview of all screwdrivers on pages 344 and 345 in the catalogue

MEG 603555



DIM ANALOG®

MEG 5.DIM.355



SLEEVE (belongs to DIM ANALOG®)

plaster model

DIM 5.XSL.A20





NT-PREFORM®

incl. coloured milling screw

MEG 10.Mi3.555.05



SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

MEG 9.S3D3.555



NT-MILL LAB ANALOG

for NT-Preform®

MEG 5.La3.555



COLOURED MILLING SCREW

for NT-Preform®

MEG 6.Sc3.555



ABUTMENT SCREW



Overview of all screwdrivers on pages 344 and 345 in the catalogue

MEG 603555



DIM ANALOG®

MEG 5.DIM.355



SLEEVE (belongs to DIM ANALOG®)
plaster model

DIM 5.XSL.A20



ABUTMENT GRIP

for NT-Holder fixture

W 12.200.000



I-BUTMENT

individual milled by NT-iDent

Not available in the USA

TZ15.TiO.Si1



MEG-SERIES

2-CONNECT®-ABUTMENT SET (STRAIGHT)

incl. 2-CONnect®-Coping and Coping-screw

GH 0,5 mm MEG-2CON 8035551S

GH 2,5 mm MEG-2CON 8035552S



2-CONNECT®-ABUTMENT

Male Part

GH 0,5 mm MEG-2CON 8035551M

GH 2,5 mm MEG-2CON 8035552M



2-CONNECT®-SCAN BODY 3D-GUIDE®

for intra + extra oral scanning

KS 61 = 1,8 mm

2CON 9.S3D1.800



2-CONNECT®-COPING

Female Part

2CON F-U



COPING-SCREW

Recommended tightening torques: 25 Ncm

KS 61



2-CONNECT®-DIM ANALOG®

KS 61 = 1,8 mm

2CON 5.DIM.180



SLEEVE (belongs to DIM ANALOG®)

plaster model

DIM 5.XSL.A30



2-CONNECT®-IMPLANT PICK-UP

integrated screw

KS 61 = 1,8 mm

2CON TR-00021.8



2-CONNECT®-LAB ANALOG

KS 61 = 1,8 mm

2CON 5.La1.800



**2-CONNECT®-ABUTMENT SET (ANGLED)**

incl. 2-CONnect®-Coping, mounting aid,
Coping-screw and Male Part Screw, angled 20°,
Gingiva Height 2,5 mm

MEG-2CON 80355A3S

**2-CONNECT®-ABUTMENT**

angled 20°, Gingiva Height 2,5 mm
Male Part incl. screw,
mounting aid

MEG-2CON 80355A3M

**2-CONNECT®-SCAN BODY 3D-GUIDE®**

for intra + extra oral scanning
KS 62 = 1,6 mm

2CON 9.S3D1.600

**2-CONNECT®-COPING**

Female Part

2CON F-U

**COPING-SCREW**

Recommended tightening torques: 15 Ncm

KS 62

**MALE PART SCREW**

for angled 20° 2-CONnect® Abutment
Recommended tightening torques: 20 Ncm

MEG 6.Ms3.555

**2-CONNECT®-DIM ANALOG®**

KS 62 = 1,6 mm

2CON 5.DIM.160



SLEEVE (belongs to DIM ANALOG®)
plaster model

DIM 5.XSL.A30

**2-CONNECT®-IMPLANT PICK-UP**

integrated screw
KS 62 = 1,6 mm

2CON TR-00031.6

**2-CONNECT®-LAB ANALOG**

KS 62 = 1,6 mm

2CON 5.La1.600

