


NT



TITANIUM BASE NX

*Titanium Base NX (NoHex) for implant
supported hybrid bridges/-bars*

NT-TRADING TIBASE NX LIBRARY

Non-engaged titanium bases (NoHex) for the production of prosthetic constructions on multiple implants (bars, bridges). These are available in 2 different gingival heights (GH) and can be customized in height from 5mm (CH5) to 3mm (CH3).

BENEFITS:

- non-engaged
- cylinder height can be adjusted from 5 mm to 3 mm
- different gingiva heights available
- perfect adaption to the gingival margin by optimized transition zone between implant and abutment
- improved for CAD/CAM process

APPLICATION:

- for individual hybrid bridges/ -bars on multi-implant cases made from ceramics (zirconium oxide, lithium disilicate) or synthetic (PMMA, PEEK, composite)



TECHNICAL DETAILS:

- Titanium Grade 5 ELI




TITANIUM BASE NX

SELECTION OF TIBASE NX (NOHEX):

3shape



exocad



ITEM SELECTION (FAX TEMPLATE)


SERIES	PLATFORM	1. GINGIVA HEIGHT	ITEM NUMBER	2. GINGIVA HEIGHT	ITEM NUMBER	QUANTITY 1. GH 2. GH	
BEG-SERIE  Soon available	Ø 3.25/3.75	GH 0.3 mm	BEG82323703	GH 1.3 mm	BEG82323713		
	Ø 4.1	GH 0.3 mm	BEG82410003	GH 1.3 mm	BEG82410013		
	Ø 4.5	GH 0.3 mm	BEG82450003	GH 1.3 mm	BEG82450013		
	Ø 5.5	GH 0.2 mm	BEG82550002	GH 1.2 mm	BEG82550012		
BIO-SERIES  	Ø 3.0	GH 0.5 mm	BIO 82300005	GH 1.5 mm	BIO 82300015		
	Ø 3.5	GH 0.5 mm	BIO 82350005	GH 1.5 mm	BIO 82350015		
	Ø 4.5	GH 0.5 mm	BIO 82450005	GH 1.5 mm	BIO 82450015		
	Ø 5.7	GH 0.5 mm	BIO 82570005	GH 1.5 mm	BIO 82570015		
BRE-SERIES  Soon available	Ø 3.5	GH 0.5 mm	BRE 82350005	GH 1.5 mm	BRE 82350015		
	Ø 4.0	GH 0.5 mm	BRE 82400005	GH 1.5 mm	BRE 82400015		
CAM-SERIES  	Ø 3.3	GH 0.5 mm	CAM 82330005	GH 1.5 mm	CAM 82330015		
	Ø 3.8	GH 0.5 mm	CAM 82380005	GH 1.5 mm	CAM 82380015		
	Ø 4.3	GH 0.5 mm	CAM 82430005	GH 1.5 mm	CAM 82430015		
	Ø 5.0	GH 0.5 mm	CAM 82500005	GH 1.5 mm	CAM 82500015		
COL-SERIES  	Ø 3.3	GH 1.0 mm	COL 82330010	GH 2.0 mm	COL 82330020		
	Ø 3.8/4.3	GH 1.0 mm	COL 82384310	GH 2.0 mm	COL 82384320		
	Ø 5.0	GH 1.0 mm	COL 82500010	GH 2.0 mm	COL 82500020		
E-SERIES  	Ø 3.5 NP	GH 0.4 mm	E 82350004	GH 1.4 mm	E 82350014		
	Ø 4.3 RP	GH 0.4 mm	E 82430004	GH 1.4 mm	E 82430014		
	Ø 5.0 WP	GH 0.4 mm	E 82500004	GH 1.4 mm	E 82500014		
	Ø 6.0	GH 0.4 mm	E 82600004	GH 1.4 mm	E 82600014		

fax form


SERIES	PLATFORM	1. GINGIVA HEIGHT	ITEM NUMBER	2. GINGIVA HEIGHT	ITEM NUMBER	QUANTITY 1. GH 2. GH	
F-SERIES  	Ø 3.5 NP	GH 0.8 mm	F 82350008	GH 1.8 mm	F 82350018		
	Ø 4.3/5.0 RP	GH 0.8 mm	F 82435008	GH 1.8 mm	F 82435018		
H-SERIES  	Ø 3.4	GH 0.3 mm	H 82340003	GH 1.3 mm	H 82340013		
	Ø 4.1	GH 0.3 mm	H 82410003	GH 1.3 mm	H 82410013		
	Ø 5.0	GH 0.3 mm	H 82500003	GH 1.3 mm	H 82500013		
HIO-SERIES  Soon available	Ø 3.5 MC	GH 1.0 mm	HIO 82350010	GH 2.0 mm	HIO 82350020		
	Ø 4.0-6.0 RC	GH 1.0 mm	HIO 82406010	GH 2.0 mm	HIO 82406020		
K-SERIES  Soon available	Ø 3.5 NP	GH 0.6 mm	K 82350006	GH 1.6 mm	K 82350016		
	Ø 4.1 RP	GH 0.5 mm	K 82410005	GH 1.5 mm	K 82410015		
	Ø 5.1 WP	GH 0.6 mm	K 82510006	GH 1.6 mm	K 82510016		
L-SERIES  	Ø 3.3 NC	GH 0.7 mm	L 82330007	GH 1.7 mm	L 82330017		
	Ø 4.1/4.8 RC	GH 0.7 mm	L 82414807	GH 1.7 mm	L 82414817		
MEG-SERIES  Soon available	Ø 3.5-5.5	GH 0.3 mm	MEG 82355503	GH 1.3 mm	MEG 82355513		

fax form

ITEM SELECTION (FAX TEMPLATE)

SERIES	PLATFORM	1. GINGIVA HEIGHT	ITEM NUMBER	2. GINGIVA HEIGHT	ITEM NUMBER	QUANTITY 1. GH 2. GH	
MEX-SERIES 	Ø 3.7-4.8 STD	GH 1.0 mm	MEX 82414810	GH 2.0 mm	MEX 82414820		
MIS-SERIES  Soon available	Ø 3.3 NP Ø 3.75/4.2 SP Ø 5.0 WP	GH 0.5 mm	MIS 82330005 MIS 82374205 MIS 82500005	GH 1.5 mm	MIS 82330015 MIS 82374215 MIS 82500015		
N-SERIES  Soon available	Ø 4.8 RN Ø 6.5 WN	GH 0.2 mm GH 0.3 mm	N 82480002 N 82650003	– –	– –	–	–
						–	–
NEO-SERIES  Soon available	Ø 3.25 Ø 3.5-5.5	GH 0.5 mm GH 0.5 mm	NEO 82320005 NEO 82355505	GH 1.5 mm GH 1.5 mm	NEO 82320015 NEO 82355515		
R-SERIES 	Ø 3.5 Ø 4.5 Ø 5.7	GH 0.5 mm GH 0.5 mm GH 0.5 mm	R 82350005 R 82450005 R 82570005	GH 1.5 mm GH 1.5 mm GH 1.5 mm	R 82350015 R 82450015 R 82570015		
S-SERIES  Soon available	Ø 3.5/4.0 Ø 4.5/5.0	GH 0.8 mm GH 0.8 mm	S 82354008 S 82455008	GH 1.8 mm GH 1.8 mm	S 82354018 S 82455018		

fax form

SERIES	PLATFORM	1. GINGIVA HEIGHT	ITEM NUMBER	2. GINGIVA HEIGHT	ITEM NUMBER	QUANTITY 1. GH 2. GH	
SEV-SERIES 	Ø 3.0 Ø 3.6 Ø 4.2 Ø 4.8 Ø 5.4	GH 1.0 mm GH 1.0 mm GH 1.0 mm GH 1.0 mm GH 1.0 mm	SEV 82300010 SEV 82360010 SEV 82420010 SEV 82480010 SEV 82540010	GH 2.0 mm GH 2.0 mm GH 2.0 mm GH 2.0 mm GH 2.0 mm	SEV 82300020 SEV 82360020 SEV 82420020 SEV 82480020 SEV 82540020		
T-SERIES 	Ø 3.4 Ø 3.8 Ø 4.5 Ø 5.5	GH 0.3 mm GH 0.3 mm GH 0.6 mm GH 0.6 mm	T 82340003 T 82380003 T 82450006 T 82550006	GH 1.3 mm GH 1.3 mm GH 1.6 mm GH 1.6 mm	T 82340013 T 82380013 T 82450016 T 82550016		
Y-SERIES  Soon available	Ø 3.5 C/	GH 1.0 mm	Y 82350010	GH 2.0 mm	Y 82350020		

fax form

- BEG-Series compatible to BEGO Implant Systems®Semados®S, RS-RX, RI-Line
- BIO-Series compatible to BioHorizons ®/ internal
- BRE-Series compatible to Bredent Medical®Sky®
- CAM-Series compatible to Altatec®/CAMLOG®
- COL-Series compatible to Altatec®/Conelog®
- E-Series compatible to NOBEL BIO CARE™/REPLACE SELECT®
- F-Series compatible to NOBEL BIO CARE™/NOBEL ACTIVE™
- H-Series compatible to BIOMET 3i®/CERTAIN®
- HIO-Series compatible to HIOSSEN®OSSTEM TS®
- K-Series compatible to NOBEL BIO CARE™Brånemark®
- L-Series compatible to STRAUMANN/BONELEVEL®
- MEG-Series compatible to MEGAGEN /AnyRidge®
- MEX-Series compatible to Medentis Medical®/ ICX®
- MIS-Series compatible to MIS®SEVEN®
- NEO-Series compatible to NEOSS Implant System®NEOSS™
- N-Series compatible to Straumann / SynOcta®
- R-Series compatible to ZIMMER Dental/ TAPERED SCREW-VENT®
- S-Series compatible to DENTSPLY Implants®Astra Tech Implant SystemOsseoSpeed® TX
- SEV-Series compatible to DENTSPLY Implants®/ASTRA TECH Implant System™ EV
- T-Series compatible to DENTSPLY Implants®/XIVE®
- Y-Series compatible to DENTSPLY Implants® /ANKYLOS®



nt-trading GmbH & Co. KG

G.-Braun-Straße 18 · 76187 Karlsruhe · Germany

info@nt-trading.com · www.nt.dental
phone: +49 (0) 721-915471-60 · fax: +49 (0) 721-915471-61